

SPECIFICATION

Mycobacterium tuberculosis MPT64 (TB-MPT64) Antibody

GENERAL INFORMATION

Catalogue Number	HM1423	Storage Temperature	Short Term: 2-8°C. Long Term: -20°C.
		Shipping Note	Cold Packs

PRODUCT INFORMATION

Format	Purified Monoclonal Antibody
Host/Source	Mouse
Purity	>95% by SDS-PAGE
Immunogen/Antigen	M. tuberculosis MPT64 (recombinant)
Class/Isotope	IgG
Supplied In	Phosphate Buffered Saline, pH 7.0
Preservative	0.02% Sodium Azide
Applications	Lateral Flow
Matched Pair	HM1424 (detection)

HANDLING & SAFETY INFORMATION

Handling	Avoid repeated freezing and thawing.
----------	--------------------------------------

