

SPECIFICATION

Rheumatoid Factor Concentrate Human Plasma

GENERAL INFORMATION

Catalogue Number	508-05	Storage Temperature	-20°C
Synonyms	RF, RF Concentrate	Storage Note	Avoid multiple freeze/thaws
		Shipping Note	Dry Ice

PRODUCT INFORMATION

Source	Human Defibrinated Plasma (multiple, pooled donors)
Form	Frozen Solution
Appearance	Yellow-Orange-Brown Liquid
Activity	≥ 1500 IU/mL (Rheumatoid Factor Turbidimetric Immunoassay)
Protein	Report in mg/mL (Biuret)
Specific Activity	Report in IU/mg protein
Recertification	One year from date of manufacture

HANDLING INFORMATION

Handling	Handle as a potentially hazardous substance.
Recovery	If small volume, microfuge vial before drawing from container.

Non-infectious disease statement: Negative or non-reactive at the donor level for anti HIV 1 and 2, anti HCV, HBsAg, HCV NAT, HIV-1 NAT, and syphilis by FDA approved methods. However, because no test method can offer complete assurance that infectious agents are absent, this material should be handled at the Biosafety Level 2 (BSL 2) as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Biosafety in Microbiological and Biomedical Laboratories", 2009. This product is intended for research, laboratory use, or further manufacturing. It is not for administration to humans. We certify this to be true to the best of our knowledge. The statements herein are offered for informational purposes only and are intended to be used solely for your consideration, investigation and verification. The suitability of this product for each application must be determined by each end user.

